



GE Healthcare



**Right when
it matters.
Right where
it counts.**



Venue Go

Venue Go™ joins the Venue™ family,* the first AI-enabled family of Point of Care ultrasound systems. It features a uniquely adaptable design that goes from cart to table to wall and shares a common platform with Venue, offering the same user interface, interchangeable probes,** strict cyber security measures, and consistent software updates to help ensure user familiarity with the systems.

gehealthcare.com

*Venue family, as referenced herein, includes Venue and Venue Go systems
**Applies to Venue R2.5 and Venue Go

Inspired by user needs

Designed for Point of Care.

Up to two-hour battery

Easy-to-see countdown timer changes color when you are running low

15-inch monitor

Sized for the environment and designed for easy repositioning and simple transport

Button probe

Allows you to directly control the system and simultaneously hold a needle

Seamless flat display

Enables thorough cleaning to support your infection-control efforts

Uniquely adaptable

Exclusive cradle enables system to go from cart to table to wall



Auto VTI

Ready where you are

Customizable LED indicators help ensure availability status

Multipurpose handle

Adjustable system height and additional cord management channels

Designed for cable management

Transducers on top keep cables off the ground

Cardiac assessment

TEE and ECG functionality

Full-screen mode

Simple user interface


Backed by a multiyear warranty

Support you can count on

Powered by AI

Auto tools for many applications.

Simplify and accelerate patient assessments with AI-enabled auto tools, available exclusively on Venue family systems. Utilizing proprietary algorithms, we synthesized the data from numerous patients to ensure accurate calculations for clinical confidence.



Auto VTI


Calculate VTI and CO in one simple step. The VTI Trending function helps clinicians quickly visualize the trend so the next course of action can be determined.



Experience 82% time savings¹



Study found better correlation than with manual CO measurements²




eFAST tool

Assess and document patient status by mapping key areas of the body.



Reduce keystrokes by 80% compared to a traditional eFAST exam³




Auto B-line tool

In one single step, this tool calculates the overall lung score by highlighting and counting B-lines in real time and displaying the image with the highest B-line count.



Study found the tool to be **comparable and as highly reliable as visual counting** performed by experts⁴



Auto IVC

Measure IVC collapsibility or distensibility accurately and automatically.



Study found the IVC measures were **equivalent to an expert user's ability 90% of the time** for minimal diameters and 97% for maximal diameters⁵

Learn more at gehealthcare.com/VenueGo



GE Healthcare

Disclaimers and references

Auto VTI – Calculate VTI and CO in one simple step. After running the Auto VTI tool, VTI Trending helps clinicians quickly visualize the trend and determine a next course of action in treatment. Auto VTI can provide up to 90% reduction in keystrokes and take up to 82% less time than manual method of calculations, as performed by experts.¹

A recent study determined in an experimental model of hemorrhagic shock by Bobbia, et al., Venue Auto VTI tool was found to be better correlated with CO measured by thermodilution than manual echocardiographic measurements.²

- 1 Based on a GE internal study with Venue GO DOC2254811.
- 2 Xavier Bobbia; Laurent Muller, et al. A New Echocardiographic Tool for Cardiac Output Evaluation: An Experimental Study 2018 OI: 10.1097/SHK.0000000000001273, PMID: 30300317

eFAST tool – Requiring up to 80% fewer keystrokes,³ this tool helps clinicians quickly assess and document patient status—from internal bleeding to a pneumothorax—by mapping key areas of the body in a way that is intuitive and aligned with clinical workflows.

- 3 eFAST Comparison Study: Manual vs. Venue Automation. GE internal study. (DOC2222911)

Auto B-line tool – This tool highlights and counts B-lines in real time, and automatically displays the image with the highest B-line count. A recent study found the Auto B-line tool to be comparable to and as highly reliable as visual counting performed by experts.⁴

- 4 Short J, Acebes C, Rodriguez-de-Lema G, et al. Visual versus automatic ultrasound scoring of lung B-Lines: reliability and consistency between systems. Med Ultrasonography 2019, Vol. 21 no. 1, 45-49 DOI: 10.11152/mu-1885

Auto IVC – Provides the clinician with the ability to measure IVC collapsibility or distensibility accurately and automatically. In one study, the IVC measures were equivalent to an expert user's ability 90% of the time for minimal diameters and 97% for maximal diameters.⁵

- 5 Venue Go R2 Technical Product Claims Document (DOC2199650)

Imagination at work

©2019 General Electric Company

All rights reserved. GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Venue Go has been tested for use in conditions similar to indoor environments and has not be evaluated for use during transport or outdoors. Use is restricted to environmental properties described in the user manual. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, Venue, Venue Go, and InSite are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB66516XX