

Codonics Grayscale Paper

Reduce hardcopy costs without sacrificing image quality using Codonics exclusive DirectVista® Grayscale Paper



“Paper printing is a boon to radiology, not only in reducing media price but also by realizing significant savings in workflow and other soft-costs.”

John C. Litchney, MBA, RT, Administrator for The Cleveland Clinic, Department of Regional Radiology

“This paper-like media is significantly more affordable than film and offers many advantages to the end user.”

Dr. David Stemerman, Radiologist, Open High-Field MRI and CT of Westchester, NY

11" x 14"
28 x 36 cm

14" x 17"
35 x 43 cm

A and
A4 Size



- Superior quality prints at a fraction of the cost of film
- Image quality is far superior to laser and office printers
- Ideal for referring physicians and patient consultations
- Letter size allows it to fit conveniently into a referring physician report or patient's chart for quick reference
- No toner, wax or ribbons to ever replace




Codonics DirectVista® Grayscale Paper



DirectVista® Grayscale Paper is specifically created for medical applications. This high quality media is ideal for referring physicians who prefer to view images in room light without a light box.

Exclusively for Codonics imagers, DirectVista Paper is available in a variety of sizes, including 14" x 17" and 11" x 14", providing cost-saving alternatives to large format film. Our A and A4 size output fits conveniently into a chart so that the image can follow the patient and doctors can easily flip to it, like a note or report.

DirectVista Paper produces stunning grayscale prints with outstanding contrast ratios. There is no toner, wax or ribbon to replace – ever. It's the ideal media for diagnostic review, patient and referral copies as well as archiving.



Reduce your film budget with DirectVista Paper. Substituting 14"x17" paper for large format film immediately reduces your hardcopy cost by 25%, while placing key images on letter size paper increases the savings up to 65% over film!



DirectVista Grayscale Paper has a heavy, photo-weight base material that resists water and tearing. Since there is no ink or toner, images never smear or transfer onto hands or other prints.

Specifications:

- Compatible with Horizon and NP/EP multi-media imagers
- Continuous tone grayscale rendering
- Dmin < 0.10 OD, Dmax > 2.00 OD (reflective)
- Thickness ~ 0.2mm


Available Sizes:

Horizon Imagers

- 14" x 17" (35 x 43 cm)
- 11" x 14" (28 x 35 cm)
- A (8.5" x 11")
- A4 (210mm x 297mm)

NP/EP Imagers

- A (8.5" x 11")
- A4 (210mm x 297mm)



For more information contact Codonics at **800-444-1198** or **440-243-1198** and ask for our Media Sales Department.



17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com